

# MewRecorder: Simplified Code Structure and Logic Summary (with OptiTrack)

## Overview

**MewRecorder.mlapp** is a MATLAB App Designer GUI for synchronized **ultrasound (Telemed EchoWave II)**, **audio**, and now **OptiTrack motion tracking** recording. It captures and synchronizes multimodal data, then optionally converts ultrasound `.tvd` files into combined `.mp4` video using **FFmpeg**.

---

## Main Functions

### 1. Initialization ( `startupFcn` → `initialize` )

- Sets up UI and verifies MATLAB version.
- Detects and connects to EchoWave II (Telemed ultrasound).
- Initializes audio device.
- Checks FFMPEG installation.
- Prepares optional OptiTrack recording control (if enabled).
- Sets app status: `inactive`, `ready`, or `busy`.

**Key calls:** - `DetectEchoWave()` → Connect to EchoWave. - `selectAudioDevice()` → Configure audio input. - `CheckIfAudioDeviceSettingSupported()` → Validate sampling rate and channel count.

---

### 2. Recording Workflow

#### (a) Start Recording ( `StartRecordButtonPushed` )

1. Checks for active EchoWave and probe (if ultrasound selected).
2. Enables optional **OptiTrack recording** (future logic can start a NatNet or Motive stream here).
3. Sets state to `busy` and generates timestamped filename via `GenFilename()`.
4. Starts **audio recorder** using `audiorecorder`.
5. Starts **Telemed recording** using EchoWave commands.
6. Continuously updates UI (record duration, lamp signals) until:
7. Stop pressed, OR
8. Duration/timeout reached.

#### (b) Stop Recording ( `Bn_StopPushed` )

- Stops audio, EchoWave, and (if applicable) OptiTrack recording.
  - Saves `.wav` (normalized audio), `.tvd` (ultrasound cine), and optionally OptiTrack data file.
  - Updates `_MewRecorder.log` with timestamps, filename, and annotation.
-

### 3. Data Conversion

#### Menu → Tools → Convert TVD to MP4 (via AVI)

- Converts all recorded `.tvd` ultrasound files using EchoWave CLI.
- Extracts synchronized audio from `.wav` via sync pulses.
- Merges audio and video with FFmpeg if installed.

**Core functions:** - `ConvertTVDtoVideo()` → Convert `.tvd` to AVI/MP4. - `TrimTelemedAudio()` → Aligns Telemed sync channel. - `ffmpegCombineVideoAudio()` → Combines audio/video.

---

### 4. OptiTrack Integration

**New Component:** `Ck_ifRecordOptitrack` (checkbox)

- Enables or disables OptiTrack motion tracking recording.
  - Intended to synchronize with audio and ultrasound capture.
  - Future versions may call external MATLAB or Python scripts to record via **NatNet SDK** or **Motive API**.
  - Placeholder hook exists in the recording logic to extend later.
- 

### 5. Utility Functions

Function	Description
<code>appSetStatus</code>	Updates GUI state.
<code>do_log</code>	Appends session info (filename, times, notes).
<code>CheckOutputFolderAccess</code>	Validates output folder permissions.
<code>FindTelemedEchoWavePath</code>	Locates EchoWave installation.
<code>hasFFmpegInstalled</code>	Checks if FFmpeg is in PATH.
<code>NextAvailableFileName</code>	Avoids overwriting files.

---

### 6. Key UI Components

Element	Purpose
<b>Record Telemed Ultrasound</b>	Enable ultrasound video acquisition.
<b>Record OptiTrack</b>	Enable OptiTrack 3D motion tracking (new).
<b>Normalize Audio</b>	Normalize waveform before saving.

Element	Purpose
<b>Start / Stop Buttons</b>	Begin and end multimodal recording.
<b>Audio Lamps (1–4)</b>	Show per-channel audio activity.
<b>Filename Controls</b>	Auto-generate timestamped filenames.
<b>Tools / Info Menus</b>	Conversion, preset loading, and info display.

## File Outputs

Type	Extension	Description
Audio	<code>.wav</code>	Normalized audio signal.
Ultrasound	<code>.tvd</code>	Telemed cine file.
OptiTrack	<code>.csv</code> / <code>.mat</code> ( <i>planned</i> )	3D marker trajectories.
Log	<code>_MewRecorder.log</code>	Session metadata.
Video	<code>.mp4</code>	Combined ultrasound + audio.

## Dependencies

- MATLAB R2022b or newer.
- Windows (Admin mode recommended).
- **Telemed EchoWave II** with `AutoInt1Client.dll`.
- **FFmpeg** installed in system PATH.
- **OptiTrack NatNet SDK** (*for future full integration*).

## Simplified Workflow

```

Start App → Initialize (Audio, EchoWave, OptiTrack)
↓
Press "Record" → Start audio, ultrasound, (and OptiTrack)
↓
Stop/AutoStop → Save .wav, .tvd, .csv/.mat
↓
Convert → Merge to MP4 (optional)
↓
Log saved in _MewRecorder.log

```

---

## Notes

- The OptiTrack checkbox adds multimodal flexibility for 3D motion + ultrasound research.
- Sync can be managed by time alignment or external trigger channels.
- The GUI design remains modular for easy future integration of new sensors.